



# BMM WESTON

Sterilizers



## USER MANUAL POROUS LOAD Series 1/2/3

**AB3066 Rev. 2 20/04/2023**  
Original Instructions

BMM WESTON LTD, WESTON WORKS, FAVERSHAM, KENT, ENGLAND

Head Office: Phone 01795 533441 Fax 01795 538891 Email: [mail@bmmweston.com](mailto:mail@bmmweston.com) Website: [www.bmmweston.com](http://www.bmmweston.com)  
Service Dept: Phone 01795 597127 Fax 01795 538891 Email: [service@bmmweston.co.uk](mailto:service@bmmweston.co.uk)



## Contents

|  |    |
|--|----|
| Preface .....  | 4  |
| Intended Use and Sterilization Cycles .....                    | 5  |
| (a) Porous Load 134°C Cycle .....                              | 5  |
| (b) Porous Load 121°C Cycle .....                              | 5  |
| (c) Bowie-Dick Test .....                                      | 5  |
| (d) Leak Rate Test.....  | 6  |
| (e) Air Detector Function Test.....                            | 6  |
| Porous Load Cycle.....   | 7  |
| Safety Instructions .....                                      | 9  |
| Safety Nomenclature .....                                      | 9  |
| Operating the Sterilizer .....                                 | 9  |
| Loading and Unloading the Sterilizer .....                     | 10 |
| Processed Items and Cycles.....                                | 10 |
| Door Operation .....   | 10 |
| End of Cycle.....  | 10 |
| Environmental Considerations .....                             | 10 |
| Maintenance and Management.....                                | 11 |
| Disposal Instructions.....                                     | 11 |
| Operator Controls .....  | 12 |
| Touchscreen keypad .....                                       | 13 |
| Daily Maintenance .....  | 15 |
| Door Seal.....   | 15 |
| Cleaning.....  | 15 |
| Daily Tests.....   | 15 |
| Loading Equipment.....   | 16 |
| Carriage Loading .....   | 17 |
| Loading Shelves .....  | 18 |
| Operation and Resetting Emergency Stop Button .....            | 19 |
| Starting a Cycle .....   | 20 |
| Ending a Cycle .....   | 25 |
| Auxiliary Menus: .....   | 26 |
| Supervisor Menu – Cycle Dependent Options:.....                | 27 |
| Supervisor Menus – Copying Profiles From Cycle to Cycle: ..... | 28 |
| Supervisor Menu - Menu Entry Logs: .....                       | 29 |
| Supervisor Menu - Change I.D.s and Passcodes:.....             | 30 |
| Cycle Log Store Menu:.....                                     | 30 |
| Integral Boiler Control.....                                   | 32 |
| Printer.....   | 33 |
| Spare Parts .....  | 34 |


## Preface

THIS MANUAL CONTAINS IMPORTANT INFORMATION FOR THE OPERATOR AND SUPERVISOR IN THE OPERATION AND DAILY MAINTENANCE OF THE STERILIZER.

**WARNING - ALL OPERATORS OF THE STERILIZER SHOULD BE GIVEN APPROPRIATE TRAINING IN THE USE OF THE STERILIZER. THE OPERATOR SHOULD BE AWARE OF THE TYPE OF LOADS THE STERILIZER IS DESIGNED TO PROCESS AND THOSE THAT SHOULD NOT BE PROCESSED. ALL OPERATORS SHOULD BECOME ACQUAINTED WITH THE SECTION ENTITLED SAFETY INSTRUCTIONS CONTAINED WITHIN THIS MANUAL.**



**Warning Symbol** – Black exclamation mark symbol on a yellow background contained within a black triangular band. This symbol is included on the sterilizer label indicating that accompanying documents must be consulted before operating the sterilizer in order to find out the nature of the potential hazards and any actions which have to be taken to avoid them. Particular attention should be paid to the safety instructions section and the caution and warning statements contained throughout the manual.



**Electric Flash Hazard Symbol** – Black flash symbol on a yellow background contained within a black triangular band. Sign attached to the covers and lids of all electrical enclosures, junction boxes and operator control panel. Indicates a warning to personnel that an electrical hazard will be present on the removal of the cover, lid or panel. Only qualified maintenance and test engineers should be allowed to remove these covers, lids and panels.

## Intended Use and Sterilization Cycles

**WARNING – READ THE FOLLOWING INSTRUCTIONS CAREFULLY BEFORE ANY ATTEMPT IS MADE AT OPERATING THE STERILIZER OR FIRST PLACING THE STERILIZER INTO USE.**

**WARNING – All Operators of the sterilizer should be given appropriate training in use of the sterilizer. If coded access to the sterilizer is used, each operator should be given a code to access only those functions and operations they are required to use.**

The cycles fitted to the sterilizer will include some, if not all, of the following:

### (a) Porous Load 134°C Cycle

Uses high-temperature steam to process porous items such as linen, towels, gowns and dressings; medical and surgical equipment, instruments and utensils packaged or wrapped in porous material such as fabrics or paper. Sterilization is by means of direct contact with steam at a temperature within the range 134°C to 137°C, for a minimum holding time of 3 minutes.

As air is trapped in porous loads, an efficient and reliable air removal stage is essential prior to the sterilizing stage.

Cycles may be fitted that provide a longer drying stage time. This extended drying may be required for heavier, more difficult loads where more steam condensate/moisture has been produced or retained.

This type of process should not be used for heat-sensitive materials.

### (b) Porous Load 121°C Cycle

This cycle is identical to the 134°C cycle except that the sterilization temperature is lower at 121°C. This cycle may be more suited to loads that cannot be heated to 134°C but will withstand 121°C.

Similarly, standard and extended drying cycles may be fitted.

### (c) Bowie-Dick Test

The Bowie-Dick Test cycle **must not** be used for processing loads.

The cycle is used to check the efficacy of the sterilization cycle using specially prepared test packs. The pack will contain a chemically coated paper. The chemical will change colour if subjected to saturated steam within the range 134°C to 137°C for a maximum holding time between 3 to 3.5 minutes. If air, or non-condensable gasses are present, the colour may not change or any change may be uneven. The test therefore not only confirms time and temperature conditions have been reached, but that satisfactory air removal has taken place allowing for the even penetration of steam into the pack.

The sterilizer should not be operated in the event of a Bowie-Dick Test failure.

### (d) Leak Rate Test

The Leak Rate Test **must not** be used for processing loads.

The cycle is designed to check the air tightness of the chamber, process connections, valves and door sealing. Air leaking into a chamber during a cycle may seriously affect sterilizing efficacy, preventing saturated steam at the correct temperature reaching all surfaces of a load.

The door is closed and sealed, and the cycle started. When the pressure in the chamber is at or below 60mbarA, all valves are closed and the pressure monitored. After an initial stabilization period of 5 minutes, the pressure rise is measured over the next 10 minutes. If during this 10 minute period, the pressure rises by more than 13mbar, the test will fail.

This test would normally be performed by a maintenance engineer.

The sterilizer should not be operated in the event of a Leak Rate Test failure.

### (e) Air Detector Function Test

The Air Detector Function Test **must not** be used for processing loads.

This cycle checks that the air detector is working correctly. The cycle is identical to a 134°C porous load cycle but during the air removal stage, air is introduced into the chamber at a controlled rate. This controlled rate of air leakage would not seriously affect sterilizing efficacy but shall be detected by the air detector and fail the cycle. This provides an assurance that during a normal cycle, the air detector will detect air and fail the cycle before the sterilizing efficacy is affected.

NOTE - FOR ADDITIONAL INFORMATION ON PROCESS AND TEST CYCLES REFER TO THE DEPARTMENT OF HEALTH PUBLICATION – CFPP 01-01 Parts B and C.

CFPP part 3 also provides a list that details the frequency the test cycles should be carried out:

#### **Daily test – User**

*Bowie-Dick test for steam penetration*


#### **Weekly tests – CP(D)**

1. *Weekly safety checks*
2. *Air leakage test*
3. *Air detector function test*
4. *Automatic control test*
5. *Bowie-Dick test for steam penetration*

## Porous Load Cycle

The pressure profile for the stages of a process cycle is shown below and consists of:

- Air Removal
- Sterilization
- Drying
- Venting



The Cycle is designed to remove air within sterilizer chamber and that in contact with the load, before subjecting the load to steam at the sterilizing temperature.

Air removal is achieved by means of a vacuum/steam fractionation process. Air is first removed from the chamber by applying a vacuum by means of a vacuum pump. Extraction stops, and steam is admitted to the chamber. This steam mixes with the residual air left after the initial vacuum. Extraction resumes, removing the steam and entrained air. The steam/extraction process is further repeated to effectively remove all the air. Air removal may be accomplished using sub-atmospheric or super-atmospheric pulses, or a combination of the two. When steam is finally admitted during the sterilizing stage of the cycle, the steam will now contact all surfaces of load in the absence of air. Air within or surrounding a load item would present a barrier to the steam and it likely that effective sterilization conditions would not be achieved.

Within a healthcare environment i.e. a CSSD, the goods to be sterilized will invariably be placed in procedure trays and the whole wrapped in special paper or fabric. As a consequence air removal is paramount.

Following air removal, steam is introduced into the chamber to raise the temperature of the load. The temperature is controlled within a defined sterilization temperature band. The load is held at this temperature for a defined period of time. BS EN 285:2015 clause 8.2.1.2.4 states that: *'The holding time shall be not less than 15 min and 3 min for sterilization temperatures of 121 °C and 134 °C respectively.'*

*NOTE - Other combinations of temperature and time may be required.*

## INTENDED USE- POROUS LOADS

BS EN 285:2015 clause 8.2.1.1 states that: *'The sterilization temperature band shall have the lower limit defined by the sterilization temperature and an upper limit of + 3 K. This upper limit value should not be exceeded throughout the whole operating cycle.'*

The temperature reference thermocouple is located in the active discharge line from the chamber. This vent is normally the drain line at the base of the chamber.

At the end of the sterilization stage, the load items and their wrappings will be wet. Steam will have condensed on the surfaces of the items in the raising their temperature. To maintain the sterility of wrapped items on removal from the sterilizer the wrapping must be dry. Any dampness in the wrapping could promote the growth of bacteria and possible contamination of the load. To dry the load, a vacuum is applied to the sterilizer chamber. As the pressure steadily decreases, and the boiling point of the water present is lowered the water will start to flash off as steam. For complete vaporisation the remaining water must gain heat, the specific enthalpy of vaporisation, from the surrounding heated load. This migration of heat can be slow and accounts for the drying stage being a lengthy part of the complete cycle. Typically, a vacuum level of 40mbarA, or lower, is attained and held. At this pressure, water will boil at around 30°C.

At the end of drying phase, the chamber pressure must be returned to atmospheric pressure before the sterilizer door can be opened and the items removed. This is achieved by admitting air into the chamber via a bacteria retentive filter. If bacteria were present in the air being introduced, these could pass through the wrapping and contaminate the sterilized items.

Dry wrapped items removed from the sterilizer will retain their sterility for a considerable period and will be suitable for storage in a clean environment before being re-used.

**For details of the parameters and tolerances critical to the sterilization process refer to the commissioning report or most recent re-validation report. These reports will also provide a record of the pressure/temperature profiles against time for the standard loads and cycles tested at that time.**



## Safety Instructions

**It is the responsibility of the User to ensure that all personnel who come into contact with the sterilizer fully understand and comply with the safety instructions contained within this manual. This applies equally to all personnel engaged in the operation, maintenance and testing of the sterilizer.**

### Safety Nomenclature

Instructions may be in the form of NOTES, CAUTIONS, and WARNINGS, and are used in the manual to draw attention to important work and safety practices.

**NOTE –** A NOTE IS USED TO DRAW ATTENTION TO AN INSTRUCTION THAT IS NOT DIRECTLY RELATED TO PERSONNEL SAFETY OR DAMAGE TO THE MACHINE.

**CAUTION -** A CAUTION IS USED TO INDICATE A HAZARDOUS SITUATION THAT MAY RESULT IN INJURY TO PERSONNEL OR DAMAGE TO THE MACHINE.

**WARNING -** A WARNING IS USED TO INDICATE A HAZARDOUS SITUATION WHICH HAS SOME PROBABILITY OF LEADING TO DEATH OR SERIOUS PERSONAL INJURY.

**WARNING – READ THE FOLLOWING INSTRUCTIONS CAREFULLY BEFORE ANY ATTEMPT IS MADE WHEN OPERATING THE STERILIZER FOR THE FIRST TIME OR WHEN FIRST PLACING THE STERILIZER INTO USE.**

The following instructions should be carried out in order that the sterilizer is operated safely, and that maintenance and validation tasks are undertaken at regular intervals to ensure that the sterilizer continues to function safely and product continues to be processed satisfactorily.

### Operating the Sterilizer

**WARNING – All Operators should be given appropriate training in use of the sterilizer. If coded access to the sterilizer is used, each operator should be given a code to access only those functions and operations they are authorized to use. Training should be ongoing and reviewed regularly. Training records should be documented.**

**WARNING – If the equipment is used in a manner not specified in this manual, the protection provided by the equipment may be impaired.**

**CAUTION -** Do not touch any part of the sterilizer chamber or door. These parts are likely to be hot and could cause burns if touched without wearing heat protective clothing and gloves.

**WARNING – Do not enter or partially enter the sterilizer chamber without first obtaining the key for the control panel auto/manual door switch and then turning the switch to the manual position. With the switch in the manual position, remove the key and keep on your person until it is safe to return the switch to the auto position and for normal operation of the sterilizer to resume.**

**CAUTION –** All sterilizer keys should be held in safe and secured storage. The issuing of keys should be authorised by the person/s responsible for their safe storage.

**CAUTION –** After use, no key should be left in any key switch. The key should be returned to the person/s responsible for the safe storage of the sterilizer keys.

**WARNING – An electrical power disconnection device (isolator) is sited in the plantroom area. This should be turned off and locked if work of a hazardous nature is to be performed on the sterilizer. Equipment should not be positioned in such a way that it is difficult to access this device.**

## Loading and Unloading the Sterilizer

CAUTION - Always wear heat protective gloves when handling, removing, or partially removing load carriers or shelves from the sterilizer chamber. Load carriers, shelves and load items are likely to be hot and could cause burns if touched without wearing heat protective gloves.

CAUTION – When transferring a load carrier from the sterilizer chamber onto a trolley, or from the trolley into the chamber, always ensure that the trolley locking device is engaged to the sterilizer chamber locking plate.

CAUTION – Do not overload loading carriage. This may give rise to excessive effort being required to maneuver the trolley and in transferring the loading carriage on and off the trolley.

## Processed Items and Cycles

CAUTION – Articles to be sterilized should be presented to the sterilizer in a clean state in order to reduce bioburden levels. For guidance on bioburden, refer to EU standards BS EN 17665 and EN556.

CAUTION – Heat sensitive items unable to withstand the sterilizing temperatures of cycles fitted to the sterilizer, should not be processed in the sterilizer. If in any doubt as to the suitability of the sterilizing cycle, seek advice from the person responsible for the sterilizer. Advice from the manufacturer of the item to be sterilized may need to be obtained.

CAUTION – Fluids, in sealed or unsealed containers, must not be processed in sterilizers and sterilizing cycles designed solely for processing medical devices. Hot fluids can cause an explosion hazard.

## Door Operation

CAUTION – When the Sterilizer door is closing keep clear to avoid entrapment.

CAUTION – Before closing the sterilizer door ensure that no load carrier, shelf or load item protrudes from the chamber. Failure to do so may result in damage to the protruding article.

CAUTION – When opening or closing the sterilizer door, ensure that the loading trolley or shelf is clear of the moving door. Failure to do so may result in damage to the sterilizer panelwork.

Caution – Do not allow articles to be sterilised to overhang the sides of a load carrier or shelving system. Articles or their packaging may be damaged when they are being transferred into or out of the sterilizer, and this may affect their sterile condition after removal from the sterilizer.

## End of Cycle

**WARNING – Always check the printed sterilizer cycle log at cycle end to confirm whether the cycle was a pass cycle and that the cycle was appropriate for the sterilization of the processed goods.**

**WARNING - Never release goods for use from a failed cycle – all goods should be reprocessed to achieve sterilization.**

**WARNING – If a fault condition message appears on the operator message screen, report fault condition to person/s responsible for the sterilizer so that appropriate action can be taken in determining the reason/s for the fault.**

## Environmental Considerations

CAUTION – This sterilizer has been tested and conforms to the requirements of EN 61000 standard on Electromagnetic Compatibility. If the sterilizer is sited in an area where the limits of electromagnetic radiation, applied by this standard, are likely to be exceeded, i.e. close to X-ray or magnetic resonance imaging equipment, advice should be obtained from the manufacturers of both the sterilizer and the high electromagnetic source equipment.

CAUTION – Peripheral equipment connected to the sterilizer external ports i.e. data collection computers, must conform to BS EN 60950 – Specification for safety of information technology equipment, including electrical business equipment.

CAUTION – The Sterilizer should be operated in a local environment where the temperature does not exceed 40°C and a relative humidity not exceeding 80%.

CAUTION – For the Operator screens to be read and operated correctly, the screen should not be subject to direct sunlight or other bright light sources.

CAUTION – Adequate lighting should be provided in the plantroom and the operator work area to enable indicators and gauges to be readable (by normal or corrected vision) from a distance of 1 m at any external illumination level in the range of (215 ± 15) lx to (1 500 ± 15) lx.

## Maintenance and Management

**WARNING** – Before any attempt is made to work on the sterilizer, ensure that all services to the machine are isolated and where appropriate locked, and that any pressure retaining parts are safely reduced to atmospheric pressure.

**WARNING** – The installation of the sterilizer should only be undertaken by suitably trained and skilled personnel. The installation procedures and checks should be documented to demonstrate that the installation is safe to use and has been installed to specification.

**WARNING** – A scheme of inspection for the safety of the sterilizer pressure vessel and door should be instituted. A Competent Person for the inspection of pressure vessels should be appointed.

**WARNING** – After installation, the sterilizer should be commissioned and validated to demonstrate that the sterilizer, as installed, is safe to use and operates to specification. This is a documented procedure, and should be undertaken by suitably qualified personnel.

**WARNING** – In use, the sterilizer should be periodically maintained and tested to ensure that the sterilizer continues to function safely and that its performance continues to be within specification. This is a documented procedure, and should be undertaken by suitably qualified personnel.

**WARNING** – The safe state of the equipment must be verified after any repair.

**NOTE** – For the management of Sterilizers see the requirements detailed in Department of Health guidance note Choice Framework for local Policy and Procedures 01-01, Parts A, B and C.

## Disposal Instructions

According to the WEEE Directive 2002/96/EC, all categories of electrical and electronic equipment (EEE) covered by this directive in the ANNEX IA should be disposed of and collected separately and use the best available treatment, recovery and recycling techniques.

EEE contains hazardous substances to the human health and environment but EEE is also a valuable resource of raw materials. Therefore it is important to collect WEEE separately from other wastes.

BMM Weston sterilizers are subject to the Directive and we therefore urge you to dispose of the equipment separately from 'normal' household waste and make sure that it is properly treated, recovered or recycled to protect human health and the environment.

Sterilizer packaging (wooden crates, plastic wrapping etc.) should be disposed of in the same manner as conventional waste of the same material.



### **Applies to Electrical Control system fitted to Sterilizer**


The symbol shown here means that the product is classed as Electrical or Electronic Equipment and should not be disposed with other household or commercial waste at the end of its working life.

The Waste of Electrical and Electronic Equipment (WEEE) Directive (2002/96/EC) has been put in place to recycle products using best available recovery and recycling techniques to minimize the impact on the environment, treat any hazardous substances and avoid increasing landfills.

Contact your local authority or the nearest collection site to dispose of Electrical and Electronic Equipment waste.

For disposal of entire sterilizer (end-of-life) contact BMM Weston for recommendations.

## Operator Controls



### PANEL LOCK

Fold away pull lever allowing access to door switches, door speed controllers etc.

### WARNING –

THE PULL DOWN PANEL SHOULD ONLY BE REMOVED BY SUITABLY QUALIFIED MAINTENANCE PERSONNEL. REMOVAL OF THE PANEL EXPOSES PERSONNEL TO MOVING PARTS AND CRUSHING HAZARDS.

### OPERATOR VDU

The Operator Touchscreen and VDU (visual display unit) are used for cycle selection, data input, sterilizer configuration, and maintenance tasks. In addition the screen can display menus for Supervisor, Validation and Maintenance requirements as well as door and boiler control.

### PRINTER

The printer is situated below the door auto/manual Keyswitch and emergency stop button. The printer provides a record of process variables, stages and times throughout the sterilizing cycle – the cycle log. In addition to the cycle log, the contents of the Supervisor, Validation and Configuration menus can be printed for archive purposes.

### EMERGENCY STOP

Operation of the emergency stop switch will remove electrical power from all electric motors and control solenoid valves. If the door or sealing system is operating, this will be immediately halted. No door operations are possible until the emergency stop system has been reset. *Warning – see section on the Operation and Resetting of the Emergency Stop Button.*

### DOOR AUTO/MANUAL KEYSWITCH

Allows for the door to be switched from automatic to manual control and vice-versa.

# OPERATOR CONTROLS

## Touchscreen keypad

The Operator touch sensitive screen is used for inputting data and control of the sterilizer.

Alphanumeric Keypad for data entry:

|      |   |   |   |   |   |   |   |   |   |   |                                      |
|------|---|---|---|---|---|---|---|---|---|---|--------------------------------------|
|      |   |   |   |   |   |   |   |   |   | X |                                      |
|      |   |   |   |   |   |   |   |   |   |   |                                      |
| 1    | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 0 | ← | Back Space Key                       |
| Q    | W | E | R | T | Y | U | I | O | P | ↵ |                                      |
| ↓    | Z | X | C | V | B | N | M | : | ; | ↑ | Arrow Keys for<br>Movement of Cursor |
| &123 |   |   |   |   |   |   |   |   | ← | → |                                      |

On pressing the '&123' key the following keypad will be displayed:

|     |   |   |   |   |   |   |   |   |   |     |                                      |
|-----|---|---|---|---|---|---|---|---|---|-----|--------------------------------------|
|     |   |   |   |   |   |   |   |   |   | X   |                                      |
|     |   |   |   |   |   |   |   |   |   |     |                                      |
| 1   | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 0 | ←   | Back Space Key                       |
| _   | / | : | ; | ( | ) | € | & | @ |   | ↵   |                                      |
| #+= |   | . | , | ? | ! | ' | " |   |   | #+= | Arrow Keys for<br>Movement of Cursor |
| ABC |   |   |   |   |   |   |   |   | ← | →   |                                      |

On pressing either of the '#+= ' keys the following keypad will be displayed:

|     |   |   |   |   |   |    |   |   |   |     |                                      |
|-----|---|---|---|---|---|----|---|---|---|-----|--------------------------------------|
|     |   |   |   |   |   |    |   |   |   | X   |                                      |
|     |   |   |   |   |   |    |   |   |   |     |                                      |
| 1   | 2 | 3 | 4 | 5 | 6 | 7  | 8 | 9 | 0 | ←   | Back Space Key                       |
| [   | ] | { | } | # | % | ^  | * | + | = |     |                                      |
| -   | \ |   | ~ | < | > | \$ | £ | ¥ |   | ↵   | Arrow Keys for<br>Movement of Cursor |
| 123 |   | . | , | ? | ! | '  | " |   |   | 123 |                                      |
| ABC |   |   |   |   |   |    |   |   | ← | →   |                                      |

## OPERATOR CONTROLS

Numeric Keypad for data entry:

On selecting a field, the current value will be displayed for that field.

|     |   |   |     |
|-----|---|---|-----|
|     |   |   | X   |
| min | 9 |   | max |
| 1   |   |   | 31  |
| 7   | 8 | 9 | ←   |
| 4   | 5 | 6 | ±   |
| 1   | 2 | 3 | ↵   |
| 0   |   | . |     |

Back Space Key

Addition of +/- symbol (where appropriate)

Enter Key

**min** – the minimum value that can be entered

**max** – the maximum value that can be entered

When the new value has been set, pressing the Enter Key will update the field.

## Daily Maintenance

These tasks should be carried out by the user/operator.

**CAUTION –** BEFORE CARRYING OUT ANY OF THE TASKS DETAILED BELOW, READ THOROUGHLY THE GENERAL SAFETY INSTRUCTION CONTAINED WITHIN THIS MANUAL.

### Door Seal

Clean the surface of the door seal with a damp cloth, wetted with water, to remove any deposits.

Visually inspect the condition of the door seal checking for any cuts or damage that might affect its ability to seal when the door is in a sealed state. If any doubt exists as to the condition of the door seal, report this to the person responsible for the sterilizer.

**CAUTION –** WHEN CLEANING OR INSPECTING THE DOOR SEAL, ALWAYS WEAR HEAT PROTECTIVE GLOVES. DOOR AND CHAMBER PARTS ARE LIKELY TO BE VERY HOT.

### Cleaning

The stainless steel facia panel can be cleaned with a damp cloth, wetted with water or disinfectant. A liquid stainless steel polish can be used but check that its use does not present any microbial problems.

Clean the operator VDU and printer cover with a damp cloth, wetted with water. Do not use any solvent cleaners on these surfaces.

**WARNING –** DO NOT ENTER OR PARTIALLY ENTER THE STERILIZER CHAMBER WITHOUT FIRST OBTAINING THE KEY FOR THE CONTROL PANEL AUTO/MANUAL DOOR SWITCH AND THEN TURNING THE SWITCH TO THE MANUAL POSITION. WITH THE SWITCH IN THE MANUAL POSITION, REMOVE THE KEY AND KEEP ON YOUR PERSON UNTIL IT IS SAFE TO RETURN THE SWITCH TO THE AUTO POSITION AND FOR NORMAL OPERATION OF THE STERILIZER TO RESUME.

**CAUTION –** ALL STERILIZERS KEYS SHOULD BE HELD IN SAFE AND SECURED STORAGE. THE ISSUING OF KEYS SHOULD BE AUTHORISED BY THE PERSON/S WHO IS/ARE RESPONSIBLE FOR THE SAFE STORAGE OF THE STERILIZER KEYS.

**CAUTION –** WHEN REMOVING THE CHAMBER FILTERS AND/OR CLEANING THE CHAMBER, ALWAYS WEAR HEAT PROTECTIVE GLOVES AND CLOTHING. CHAMBER PARTS ARE LIKELY TO BE VERY HOT AND TO CAUSE BURNS.

### Instrumentation

Check the VDU to establish that all the services are turned on and that the correct readings are indicated on the VDU.

### Daily Tests

A Bowie and Dick Cycle should be run daily, preferably at the beginning of the day. The chemical indicator will demonstrate the effectiveness of the air removal process. If an even colour change is not obtained across the indicator paper, do not use the sterilizer until the cause is fixed.

A Leak Rate Test should be run on a weekly, if not daily basis, to check the integrity of the chamber/door system and associated pipework. The test is designed to check the air tightness of the system and to fail the cycle if air leakage is above a pre-defined value. If the Leak Rate Test indicates a failure, inform the person responsible for the sterilizer.

**WARNING -** IF ANY OF THE REQUIREMENTS FOR THAT TEST ARE NOT MET, DO NOT USE THE STERILIZER. REPORT THE FAILURE TO THE PERSON RESPONSIBLE FOR THE STERILIZER. THE STERILIZER SHOULD ONLY BE USED AFTER THE CAUSE OF THE FAILURE HAS BEEN ESTABLISHED AND CORRECTIVE ACTION HAS BEEN TAKEN.

Check the log book and production records. Complete as required.

### Printer

Check that the printer is loaded with sufficient paper for the day's operation.

## Loading Equipment


The general arrangement of the sterilizer loading equipment is shown in the diagram below.

The equipment comprises a loading trolley upon which sits the loading carriage. The loading carriage is fitted with wheels and is guided onto and off of the trolley by means of rails fitted to the trolley.

A similar set of rails is fitted to the floor of the chamber, and when the trolley and chamber rails are aligned, the loading carriage can be transferred into or out of the chamber.

With the carriage pulled onto the trolley and the carriage wheel contacting the carriage stop, the locking lever in the locked position will raise the carriage locking bar to secure the carriage to the trolley. In this position the trolley and carriage can be wheeled safely around the work area.


**CAUTION – ALWAYS ENSURE THAT THE CARRIAGE LOCKING BAR IS IN THE RAISED POSITION, SECURING THE CARRIAGE TO THE TROLLEY BEFORE MOVING THE TROLLEY AND CARRIAGE AROUND THE WORK AREA.**



With the sterilizer door in the open position, the loaded carriage and trolley is pushed up to the chamber. The trolley should be steered to align with the chamber. As the trolley approaches the chamber, the guides fitted to the front of the trolley rails will begin to engage the chamber rails to align the two sets of rails. When the two sets of rails are abutted, the trolley locking lever is turned to the unlock position. This will cause the carriage locking bar to be lowered whilst at the same time, the trolley locking bar to engage the chamber locking plate within the chamber. The engagement of the trolley locking bar to the chamber locking plate will prevent the trolley moving whilst the carriage is either transferred into or out of the sterilizer. For additional security, the trolley wheel locks should be operated, by pushing the wheel locking bars downward.

**CAUTION – WHEN TRANSFERRING THE LOADING CARRIAGE INTO OR OUT OF THE CHAMBER, ENSURE THAT THE TROLLEY LOCKING BAR IS SECURELY ENGAGED TO THE CHAMBER LOCKING PLATE SO THAT THE TROLLEY IS HELD IN POSITION. ALSO ENSURE THAT THE TROLLEY WHEEL LOCKS ARE OPERATED TO PREVENT MOVEMENT OF THE TROLLEY DURING THE LOADING OR UNLOADING OPERATION.**





When the carriage has been pushed fully into the chamber, the trolley locking lever is turned to the carriage lock position to release the trolley locking bar from the chamber locking plate. The trolley wheel locks can be released by raising the locking bars and the trolley can then be wheeled clear of the sterilizer.

**CAUTION –** BEFORE CLOSING THE STERILIZER DOOR ENSURE THAT THE CARRIAGE DOES NOT PROTRUDE FROM THE STERILIZER CHAMBER.

**CAUTION –** WHEN OPENING OR CLOSING THE STERILIZER DOOR, ENSURE THAT THE LOADING TROLLEY IS CLEAR OF THE DOOR. FAILURE TO DO SO MAY RESULT IN DAMAGE TO THE STERILIZER PANELWORK.

To unload the carriage from the chamber, the sequence of events as described above should be followed. When unloading the carriage from the sterilizer, or handling recently sterilized load items, always wear heat protective gloves.

**CAUTION -** ALWAYS WEAR HEAT PROTECTIVE GLOVES WHEN HANDLING/REMOVING LOAD CARRIERS OR SHELVES FROM THE STERILIZER CHAMBER. LOAD CARRIERS, SHELVES AND LOAD ITEMS ARE LIKELY TO BE HOT AND COULD CAUSE BURNS IF TOUCHED WITHOUT WEARING HEAT PROTECTIVE GLOVES.

### **Carriage Loading**

To ensure that the carriage shelves are not overloaded and that the load items are dry after a sterilization cycle, instrument trays should only be placed on the shelves in single layers. The loaded items should not overhang the sides of the shelves or carriage.

**CAUTION –** DO NOT OVERLOAD CARRIAGE AND SHELVES.

**CAUTION –** DO NOT ALLOW ARTICLES TO BE STERILISED TO OVERHANG THE SIDES OF A LOAD CARRIER OR SHELVING SYSTEM. ARTICLES OR THEIR PACKAGING MAY BE DAMAGED WHEN THEY ARE BEING TRANSFERRED INTO OR OUT OF THE STERILIZER, AND THIS MAY AFFECT THEIR STERILE CONDITION AFTER REMOVAL FROM THE STERILIZER.

## Loading Shelves

If loading shelves have been supplied with the sterilizer, they may be extended and retracted by gripping the handle and pushing or pulling as appropriate.

**CAUTION – GLOVES SHOULD BE WORN WHEN LOADING OR UNLOADING THE SHELVES.**

**CAUTION – WHEN RETRACTING THE SHELVES ENSURE THAT THEY HAVE REACHED THE STOP POINT AT THE REAR OF THE CHAMBER SO THAT THEY DO NOT IMPEDE THE MOTION OF THE STERILIZER DOOR AS IT CLOSES.**

**CAUTION – SPREAD LOADS EVENLY ACROSS THE LENGTH AND WIDTH OF THE LOADING SHELF. OVERLOADING OR IMPROPER LOADING MAY CAUSE DAMAGE OR DEFORMATION TO THE SHELF OR HARM TO THE USER.**



Shelves in closed position



Shelves in open position


# EMERGENCY STOP BUTTON

## Operation and Resetting Emergency Stop Button

When the Emergency Stop is operated, all electrical motors are isolated and all control valves become inhibited.

If the door is moving open or closed the valve will remain its position and therefore keep closing or opening.

The VDU will display




### RESETTING EMERGENCY STOP


**WARNING – MAKE SURE THAT THE REASON FOR OPERATING THE EMERGENCY STOP HAS BEEN REMOVED AND IT IS SAFE TO RETURN THE STERILIZER TO NORMAL OPERATION**

In order to reset the emergency stop an authorised key holder must use the emergency stop key to release the emergency stop. This will return the machine to a safe state before returning to the menus.

VDU if the door is in the open or intermediate position



VDU if the door is in the closed position



# STARTING A CYCLE

## Starting a Cycle

**NOTE: IF THE STERILIZER HAS AN INTEGRAL BOILER REFER TO BOILER CONTROL SECTION (TO START THE BOILER BEFORE ATTEMPTING TO RUN A CYCLE)**

**BEFORE ANY ATTEMPT IS MADE TO LOAD THE STERILIZER OR START A CYCLE, THE SAFETY INSTRUCTIONS INCLUDED AT THE FRONT OF THIS MANUAL SHOULD BE READ CAREFULLY AND OBSERVED AT ALL TIMES.**

If the door is in the closed position, press the OPEN DOOR button on the control cabinet. The door, when open, will display the main menu as below, with options to access auxiliary menus or close the door:

| Main Menu – Door Open   |        |                 |                      |                    |          |
|---|--------|-----------------|----------------------|--------------------|----------|
| WARNING CHAMBER MAY BE VERY HOT   |        |                 |                      |                    |          |
| DOOR CAN BE CLOSED  |        |                 |                      |                    |          |
| <div style="background-color: #ADD8E6; padding: 5px; margin: 5px auto; width: 150px;">AUXILIARY MENU</div> <div style="background-color: #ADD8E6; padding: 5px; margin: 5px auto; width: 150px;">CLOSE DOOR</div> |        |                 |                      |                    |          |
| CHAMBER PRESSURE  | 0959mb | BOILER PRESSURE | 4.5bar               | JACKET TEMPERATURE | 135.3°C  |
| JACKET PRESSURE   | 3160mb |                 |                      | LAST CYCLE NUMBER  | 000101   |
| DOOR 1 LOADING SIDE   |        | DATE            | Wednesday 09/05/2018 | TIME               | 10:26:28 |

Transfer the load into the sterilizer, making sure that no part of the load touches the chamber walls. Read the safety instructions section at the front of this manual for the safe handling of loads in and out of the sterilizers.

Press and hold the close door button to close the door.

When the door is fully closed the cycle selection screen will be displayed as below.

| Main Menu – Door Closed   |        |  |                      |                    |          |
|---|--------|--|----------------------|--------------------|----------|
| DOOR CAN BE OPENED  |        |  |                      |                    |          |
| <div style="background-color: #ADD8E6; padding: 5px; margin: 5px auto; width: 150px;">CYCLE 1</div>   |        | <div style="background-color: #ADD8E6; padding: 5px; margin: 5px auto; width: 150px;">CYCLE 2</div>        |                      |                    |          |
| <div style="background-color: #ADD8E6; padding: 5px; margin: 5px auto; width: 150px;">CYCLE 3</div>   |        | <div style="background-color: #ADD8E6; padding: 5px; margin: 5px auto; width: 150px;">CYCLE 4</div>        |                      |                    |          |
| <div style="background-color: #ADD8E6; padding: 5px; margin: 5px auto; width: 150px;">CYCLE 5</div>   |        | <div style="background-color: #ADD8E6; padding: 5px; margin: 5px auto; width: 150px;">CYCLE 6</div>        |                      |                    |          |
| <div style="background-color: #ADD8E6; padding: 5px; margin: 5px auto; width: 150px;">CYCLE 7</div>   |        | <div style="background-color: #ADD8E6; padding: 5px; margin: 5px auto; width: 150px;">AUXILIARY MENU</div> |                      |                    |          |
| <div style="background-color: #ADD8E6; padding: 5px; margin: 5px auto; width: 200px;">OPEN DOOR</div> |        |  |                      |                    |          |
| CHAMBER PRESSURE  | 0959mb | BOILER PRESSURE  | 4.5bar               | JACKET TEMPERATURE | 135.3°C  |
| JACKET PRESSURE   | 3160mb |  |                      | LAST CYCLE NUMBER  | 000101   |
| DOOR 1 LOADING SIDE   |        | DATE   | Wednesday 09/05/2018 | TIME               | 10:26:28 |

If sub-cycles are available the sub-cycle selection menu will be displayed on selecting a cycle:

| Sub Menu – Door Closed  |        |  |                      |                    |          |
|---|--------|--|----------------------|--------------------|----------|
| DOOR CAN BE OPENED  |        |  |                      |                    |          |
| <div style="background-color: #ADD8E6; padding: 5px; margin: 5px auto; width: 150px;">CYCLE X</div>   |        | <div style="background-color: #ADD8E6; padding: 5px; margin: 5px auto; width: 150px;">CYCLE X</div>        |                      |                    |          |
| <div style="background-color: #ADD8E6; padding: 5px; margin: 5px auto; width: 150px;">CYCLE X</div>   |        | <div style="background-color: #ADD8E6; padding: 5px; margin: 5px auto; width: 150px;">CYCLE X</div>        |                      |                    |          |
| <div style="background-color: #ADD8E6; padding: 5px; margin: 5px auto; width: 150px;">CYCLE X</div>   |        | <div style="background-color: #ADD8E6; padding: 5px; margin: 5px auto; width: 150px;">CYCLE X</div>        |                      |                    |          |
| <div style="background-color: #ADD8E6; padding: 5px; margin: 5px auto; width: 150px;">CYCLE X</div>   |        | <div style="background-color: #ADD8E6; padding: 5px; margin: 5px auto; width: 150px;">AUXILIARY MENU</div> |                      |                    |          |
| <div style="background-color: #ADD8E6; padding: 5px; margin: 5px auto; width: 150px;">OPEN DOOR</div> |        |  |                      |                    |          |
| CHAMBER PRESSURE  | 0959mb | BOILER PRESSURE  | 4.5bar               | JACKET TEMPERATURE | 135.3°C  |
| JACKET PRESSURE   | 3160mb |  |                      | LAST CYCLE NUMBER  | 000101   |
| DOOR 1 LOADING SIDE   |        | DATE   | Wednesday 09/05/2018 | TIME               | 10:26:28 |

Select the required cycle by pressing the appropriate key.

# STARTING A CYCLE

If the cycle requires an I.D. and passcode in order to start, the following screen will be displayed.

| Authorising start of cycle          |                      |                 |        |                      |               |
|-------------------------------------|----------------------|-----------------|--------|----------------------|---------------|
| ID                                  | <input type="text"/> |                 |        |                      |               |
| PASSCODE                            | <input type="text"/> |                 |        |                      |               |
| <input type="button" value="EXIT"/> |                      |                 |        |                      |               |
| CHAMBER PRESSURE                    | 0959mb               | BOILER PRESSURE | 4.5bar | JACKET TEMPERATURE   | 135.3°C       |
| JACKET PRESSURE                     | 3160mb               |                 |        | LAST CYCLE NUMBER    | 000101        |
| DOOR 1 LOADING SIDE                 |                      |                 | DATE   | Wednesday 09/05/2018 | TIME 10:26:28 |

Pressing exit will return to the cycle selection menu

On entering a valid I.D. and passcode combination, the seal door screen will be displayed as below.

| Cycle Start – Seal Door                    |        |                 |                                     |                      |               |
|--|--------|-----------------|-------------------------------------|----------------------|---------------|
| CYCLE X                                    |        |                 |                                     |                      |               |
| <input type="button" value="SEAL DOOR 1"/> |        |                 | <input type="button" value="EXIT"/> |                      |               |
| CHAMBER PRESSURE                           | 0959mb | BOILER PRESSURE | 4.5bar                              | JACKET TEMPERATURE   | 135.3°C       |
| JACKET PRESSURE                            | 3160mb |                 |                                     | LAST CYCLE NUMBER    | 000101        |
| DOOR 1 LOADING SIDE                        |        |                 | DATE                                | Wednesday 09/05/2018 | TIME 10:26:28 |

Pressing exit will return to the cycle selection menu

| Cycle Start – Select Cycle Start Option    |        |                 |                                     |                      |               |
|--|--------|-----------------|-------------------------------------|----------------------|---------------|
| CYCLE X                                    |        |                 |                                     |                      |               |
| <input type="button" value="START CYCLE"/> |        |                 | <input type="button" value="EXIT"/> |                      |               |
| CHAMBER PRESSURE                           | 0959mb | BOILER PRESSURE | 4.5bar                              | JACKET TEMPERATURE   | 135.3°C       |
| JACKET PRESSURE                            | 3160mb |                 |                                     | LAST CYCLE NUMBER    | 000101        |
| DOOR 1 LOADING SIDE                        |        |                 | DATE                                | Wednesday 09/05/2018 | TIME 10:26:28 |

Pressing exit will return to the cycle selection menu

Cycle Start

# STARTING A CYCLE

## Starting a Test Cycle

If the selected cycle is a test cycle and therefore unsuitable for processing loads, the screen will be displayed as below:

| Main Menu – Door Closed                         |        |                 |                      |                    |             |
|---|--------|-----------------|----------------------|--------------------|-------------|
|   |        |                 |                      |                    | NOT STERILE |
| CYCLE X   |        |                 |                      |                    |             |
| THIS CYCLE IS NOT SUITABLE FOR PROCESSING LOADS |        |                 |                      |                    |             |
| CONTINUE ANYWAY                                 |        |                 | EXIT                 |                    |             |
| CHAMBER PRESSURE                                | 0959mb | BOILER PRESSURE | 4.5bar               | JACKET TEMPERATURE | 135.3°C     |
| JACKET PRESSURE                                 | 3160mb |                 |                      | LAST CYCLE NUMBER  | 000101      |
| DOOR 1 LOADING SIDE                             |        | DATE            | Wednesday 09/05/2018 | TIME               | 10:26:28    |

If the cycle requires an I.D. and passcode in order to start, the following screen will be displayed.

| Authorising start of cycle |                      |                 |                      |                    |          |
|----------------------------|----------------------|-----------------|----------------------|--------------------|----------|
| ID                         | <input type="text"/> |                 |                      |                    |          |
| PASSCODE                   | <input type="text"/> |                 |                      |                    |          |
| EXIT                       |                      |                 |                      |                    |          |
| CHAMBER PRESSURE           | 0959mb               | BOILER PRESSURE | 4.5bar               | JACKET TEMPERATURE | 135.3°C  |
| JACKET PRESSURE            | 3160mb               |                 |                      | LAST CYCLE NUMBER  | 000101   |
| DOOR 1 LOADING SIDE        |                      | DATE            | Wednesday 09/05/2018 | TIME               | 10:26:28 |

On entering a valid I.D. and passcode combination, the seal door screen will be displayed as below.

| Cycle Start – Seal Door |        |                 |                      |                    |          |
|-------------------------|--------|-----------------|----------------------|--------------------|----------|
| CYCLE X                 |        |                 |                      |                    |          |
| SEAL DOOR 1             |        |                 | EXIT                 |                    |          |
| CHAMBER PRESSURE        | 0959mb | BOILER PRESSURE | 4.5bar               | JACKET TEMPERATURE | 135.3°C  |
| JACKET PRESSURE         | 3160mb |                 |                      | LAST CYCLE NUMBER  | 000101   |
| DOOR 1 LOADING SIDE     |        | DATE            | Wednesday 09/05/2018 | TIME               | 10:26:28 |

# STARTING A CYCLE

Once the door is sealed the cycle start screen will be displayed. If the cycle is configured so that a start time can be entered the set time to start cycle key will be displayed.

| Cycle Start – Select Cycle Start Option |        |                 |                      |                    |          |
|---|--------|-----------------|----------------------|--------------------|----------|
| CYCLE X                                 |        |                 |                      |                    |          |
| START CYCLE                             |        |                 | EXIT                 |                    |          |
| CHAMBER PRESSURE                        | 0959mb | BOILER PRESSURE | 4.5bar               | JACKET TEMPERATURE | 135.3°C  |
| JACKET PRESSURE                         | 3160mb |                 |                      | LAST CYCLE NUMBER  | 000101   |
| DOOR 1 LOADING SIDE                     |        | DATE            | Wednesday 09/05/2018 | TIME               | 10:26:28 |

Cycle Start

## Starting a Time-Set Cycle

If the selected cycle has the time-set option enabled, the cycle start screen will be displayed as below:

| Cycle Start – Select Cycle Start Option |        |                 |                         |                    |          |
|---|--------|-----------------|-------------------------|--------------------|----------|
| CYCLE X                                 |        |                 |                         |                    |          |
| START CYCLE                             |        |                 | SET TIME TO START CYCLE |                    |          |
| START CYCLE                             |        |                 | EXIT                    |                    |          |
| CHAMBER PRESSURE                        | 0959mb | BOILER PRESSURE | 4.5bar                  | JACKET TEMPERATURE | 135.3°C  |
| JACKET PRESSURE                         | 3160mb |                 |                         | LAST CYCLE NUMBER  | 000101   |
| DOOR 1 LOADING SIDE                     |        | DATE            | Wednesday 09/05/2018    | TIME               | 10:26:28 |

On pressing the set time to start cycle key, the user will be prompted to enter a start date and time.

| Time Set Cycle Start – Check/Set Date For Timed Start |        |                 |                      |                    |          |
|---|--------|-----------------|----------------------|--------------------|----------|
| CYCLE X   |        |                 |                      |                    |          |
| START DATE 10 / 05 / 2018                             |        |                 |                      |                    |          |
| START TIME 10 / 30                                    |        |                 |                      |                    |          |
| START CYCLE   |        |                 | EXIT                 |                    |          |
| CHAMBER PRESSURE                                      | 0959mb | BOILER PRESSURE | 4.5bar               | JACKET TEMPERATURE | 135.3°C  |
| JACKET PRESSURE                                       | 3160mb |                 |                      | LAST CYCLE NUMBER  | 000101   |
| DOOR 1 LOADING SIDE                                   |        | DATE            | Wednesday 09/05/2018 | TIME               | 10:26:28 |

# STARTING A CYCLE

On entering a valid start time the option to enter standby mode becomes available.

| Time Set Cycle Start – Check/Set Date For Timed Start |        |                 |        |                      |               |
|---|--------|-----------------|--------|----------------------|---------------|
| CYCLE X   |        |                 |        |                      |               |
| TIME SET CYCLE START SET : Thursday 10/05/2018 10:30  |        |                 |        |                      |               |
| START DATE  |        | 10              | /      | 05                   | / 2018        |
| START TIME  |        | 10              | /      | 30                   |               |
| STANDBY MODE  |        |                 |        |                      |               |
| START CYCLE   |        |                 |        | EXIT                 |               |
| CHAMBER PRESSURE                                      | 0959mb | BOILER PRESSURE | 4.5bar | JACKET TEMPERATURE   | 135.3°C       |
| JACKET PRESSURE                                       | 3160mb |                 |        | LAST CYCLE NUMBER    | 000101        |
| DOOR 1 LOADING SIDE                                   |        |                 | DATE   | Wednesday 09/05/2018 | TIME 10:26:28 |

On pressing standby mode the screen will turn blank, it can be reactivated at any point by touching it.

When the time set cycle is about to begin the following screen will be displayed.


| Cycle Start – Select Cycle Start Option              |        |                 |        |                      |               |
|--|--------|-----------------|--------|----------------------|---------------|
| <b>!! WARNING – TIME SET CYCLE ABOUT TO START !!</b> |        |                 |        |                      |               |
| CHAMBER PRESSURE                                     | 0959mb | BOILER PRESSURE | 4.5bar | JACKET TEMPERATURE   | 135.3°C       |
| JACKET PRESSURE                                      | 3160mb |                 |        | LAST CYCLE NUMBER    | 000101        |
| DOOR 1 LOADING SIDE                                  |        |                 | DATE   | Wednesday 09/05/2018 | TIME 10:26:28 |

Cycle Start




## Ending a Cycle

At the end of the cycle the screen will appear as below. The cycle name and whether the cycle was a pass or a fail will be displayed. If the cycle has failed the message 'NON-STERILE LOAD IN CHAMBER' will be displayed:




The prompt for an I.D. will only appear if the cycle has been configured to require an I.D. and passcode to end.

On entering a valid I.D. and passcode combination the option to unseal the door will be available.



Read the safety instructions section at the front of this manual for the safe operation of the door following a cycle.

Once the door is unsealed the option to open the door will be available.




Transfer the load out of the sterilizer, making sure that no part of the load touches the chamber walls. Read the safety instructions section at the front of this manual for the safe handling of loads in and out of the sterilizers.

Cycle End

## Auxiliary Menus:

The Auxiliary Menu screen can be accessed from the main menu whether the door is open or closed. The screen will be displayed as below:



### 1. SUPERVISOR MENU

Set of menus for configuring the data input requirements for starting and ending a cycle. Also general tools for handling data storage facilities.

### 2. MAINTENANCE MENU

Set of menus to assist the maintenance engineer. These include door operations, input and output status of microprocessor control signals, calibration of analogue signals, and video and printer checks.

### 3. VALIDATION MENU

A set of menus used to configure the process cycles variables. These would include variables to control the steam/vacuum pulsing during the air removal stage, sterilizing temperature limits and times, and drying vacuum levels and times.

### 4. CONFIGURATION MENU

Set of menus used to configure the basic set-up characteristics for the sterilizer. These include the specifying of the number of doors, number and type of analogue channels and the number of cycles the sterilizer will support.

### 5. ERROR HISTORY

Menu displaying the date, time and details of any errors that have occurred both in and out of cycle.

### 6. CYCLE LOG STORE MENU

Set of menus for viewing cycle logs that have been store in the sterilizer control system memory. Also details how cycle logs can be reprinted and general management features on log storage.

### 7. BOILER CONTROL MENU

Menu that describes the status of the integral boiler (if fitted) and allows for the boiler to be started or stopped.

# SUPERVISOR MENU – CYCLE DEPENDENT OPTIONS

## Supervisor Menu – Cycle Dependent Options:

The Supervisor Menu can be accessed by pressing the Supervisor Menu key from the Auxiliary Menus screen and inputting a valid I.D. and passcode combination at supervisor level or above. The screen will appear as below:

Supervisor Menu

CYCLE DEPENDENT OPTIONS      GENERAL OPTIONS

EXIT

CHAMBER PRESSURE 0959mb    BOILER PRESSURE 4.5bar    JACKET TEMPERATURE 135.3°C  
 JACKET PRESSURE 3160mb    LAST CYCLE NUMBER 000101

DOOR 1 LOADING SIDE    DATE Wednesday 09/05/2018    TIME 10:26:28

Press the Exit key to return to the Supervisor Menu

Press the Cycle Dependent Options key

Supervisor Menu - Cycle Dependent Options

CYCLE 1      CYCLE 2  
 CYCLE 3      CYCLE 4  
 CYCLE 5      CYCLE 6  
 CYCLE 7      EXIT

CHAMBER PRESSURE 0959mb    BOILER PRESSURE 4.5bar    JACKET TEMPERATURE 135.3°C  
 JACKET PRESSURE 3160mb    LAST CYCLE NUMBER 000101

DOOR 1 LOADING SIDE    DATE Wednesday 09/05/2018    TIME 10:26:28

Press the Exit Key

Select a cycle

If sub-cycles are available the sub-cycle selection menu will be displayed on selecting a cycle:

Supervisor Menu - Cycle Dependent Options

CYCLE X-1      CYCLE X-2  
 CYCLE X-3      CYCLE X-4  
 CYCLE X-5      CYCLE X-6  
 CYCLE X-7      EXIT

CHAMBER PRESSURE 0959mb    BOILER PRESSURE 4.5bar    JACKET TEMPERATURE 135.3°C  
 JACKET PRESSURE 3160mb    LAST CYCLE NUMBER 000101

DOOR 1 LOADING SIDE    DATE Wednesday 09/05/2018    TIME 10:26:28

CYCLE X-X

ID REQUIRED TO START CYCLE      Yes ▾

LEVEL OF ID REQUIRED TO START CYCLE      1    2    3    4

ID REQUIRED TO END PASS CYCLE      Yes ▾

LEVEL OF ID REQUIRED TO END PASS CYCLE      1    2    3    4

ID REQUIRED TO END FAILED CYCLE      Yes ▾

LEVEL OF ID REQUIRED TO END FAILED CYCLE      0    0    3    4

TIME SET CYCLE OPTION      Yes ▾      MINI/DETAIL LOG      Mini ▾

PRINT      SAVE PROFILE TO .STX TEXT FILE      EXIT

CHAMBER PRESSURE 0959mb    BOILER PRESSURE 4.5bar    JACKET TEMPERATURE 135.3°C  
 JACKET PRESSURE 3160mb    LAST CYCLE NUMBER 000101

DOOR 1 LOADING SIDE    DATE Wednesday 09/05/2018    TIME 10:26:28

| ITEM                                     | DESCRIPTION   |
|--|---|
| ID REQUIRED TO START CYCLE               | Dropdown box selection with Yes or No determining whether or not an ID and passcode are required to start the cycle in question.  |
| LEVEL OF ID REQUIRED TO START CYCLE      | Each field corresponds to ID level 1-4. Input the corresponding digit for the ID level to be able to start the cycle or input a 0 if the ID level is unable to start the cycle.             |
| ID REQUIRED TO END PASS CYCLE            | Dropdown box selection with Yes or No determining whether an ID and passcode are required to end the cycle in question when it passes.  |
| LEVEL OF ID REQUIRED TO END PASS CYCLE   | Each field corresponds to ID level 1-4. Input the corresponding digit for the ID level to be able to end the passing cycle or input a 0 if the ID level is unable to end the passing cycle. |
| ID REQUIRED TO END FAILED CYCLE          | Dropdown box selection with Yes or No determining whether an ID and passcode are required to end the cycle in question when it fails.   |
| LEVEL OF ID REQUIRED TO END FAILED CYCLE | Each field corresponds to ID level 1-4. Input the corresponding digit for the ID level to be able to end the failing cycle or input a 0 if the ID level is unable to end the failing cycle. |
| TIME SET CYCLE OPTION                    | Dropdown box selection with Yes or No determining whether or not the option to set a start time for the cycle in question is available.   |
| MINI/DETAIL LOG                          | Dropdown box selection with Mini or Detail determining the level of detail expressed in the data log.   |

Press the Exit key once desired changes have been made to save them. Press the Print key to print the Supervisor profile for the cycle or press the Save Profile key to save them to a text file.

# CYCLE LOG STORE MENU

## Supervisor Menus – Copying Profiles From Cycle to Cycle:

The Supervisor Menu can be accessed by pressing the Supervisor Menu key from the Auxiliary Menus screen and inputting a valid I.D. and passcode combination at supervisor level or above. The screen will appear as below:

Supervisor Menu

CYCLE DEPENDENT OPTIONS

GENERAL OPTIONS

EXIT

|                     |        |                 |                      |                    |          |
|---------------------|--------|-----------------|----------------------|--------------------|----------|
| CHAMBER PRESSURE    | 0959mb | BOILER PRESSURE | 4.5bar               | JACKET TEMPERATURE | 135.3°C  |
| JACKET PRESSURE     | 3160mb |                 |                      | LAST CYCLE NUMBER  | 000101   |
| DOOR 1 LOADING SIDE |        | DATE            | Wednesday 09/05/2018 | TIME               | 10:26:28 |

Press the Exit key to return to the Supervisor Menu

Press the General Options key

Supervisor Menu - General Options

COPY PROFILE CYCLE TO CYCLE

MENU ENTRY LOGS

CHANGE ID'S & PASSCODES

VISUAL CHART RECORDING MODE

1. ENABLE VISUAL CHART RECORDING

↓

EXIT

|                     |        |                 |                      |                    |          |
|---------------------|--------|-----------------|----------------------|--------------------|----------|
| CHAMBER PRESSURE    | 0959mb | BOILER PRESSURE | 4.5bar               | JACKET TEMPERATURE | 135.3°C  |
| JACKET PRESSURE     | 3160mb |                 |                      | LAST CYCLE NUMBER  | 000101   |
| DOOR 1 LOADING SIDE |        | DATE            | Wednesday 09/05/2018 | TIME               | 10:26:28 |

The Visual Chart Recording Mode dropdown box allows the user to enable or disable visual chart recording during the cycle by selecting the corresponding option.

Copy Supervisor Profiles Cycle to Cycle

SELECT SOURCE PROFILE (COPY FROM)

X

-

X

SELECT TARGET PROFILE (COPY TO)

X

-

X

COPY

EXIT

|                     |        |                 |                      |                    |          |
|---------------------|--------|-----------------|----------------------|--------------------|----------|
| CHAMBER PRESSURE    | 0959mb | BOILER PRESSURE | 4.5bar               | JACKET TEMPERATURE | 135.3°C  |
| JACKET PRESSURE     | 3160mb |                 |                      | LAST CYCLE NUMBER  | 000101   |
| DOOR 1 LOADING SIDE |        | DATE            | Wednesday 09/05/2018 | TIME               | 10:26:28 |

The Copy Supervisor Profile menu allows a Supervisor Profile to be copied from one cycle to another.

On pressing the Select Source Profile data entry field a keypad will appear and allow the number of the source profile to be entered in the format (x-x).

On pressing the Select Target Profile data entry field a keypad will appear and allow the number of the target profile to be entered in the format (x-x).


On pressing copy the supervisor profile data values of the source profile will be copied over to the target profile.

Press the Exit key to return to the General Options Menu.

# CYCLE LOG STORE MENU


## Supervisor Menu - Menu Entry Logs:

The Supervisor Menu can be accessed by pressing the Supervisor Menu key from the Auxiliary Menus screen and inputting a valid I.D. and passcode combination at supervisor level or above. The screen will appear as below:



Press the Exit key to return to the Supervisor Menu

Press the General Options key



Press the Menu Entry Logs Key



The screen displays a dated and time list of all supervisor menu entries and exits and their corresponding I.D.s.

Entries are coloured in green whilst exits are coloured in red.


The dropdown menu allows the corresponding information to be viewed for the Validation, Maintenance and Configuration menus.

Press the Exit key to return to the General Options Menu.


# CYCLE LOG STORE MENU

## Supervisor Menu - Change I.D.s and Passcodes:


The Supervisor Menu can be accessed by pressing the Supervisor Menu key from the Auxiliary Menus screen and inputting a valid I.D. and passcode combination at supervisor level or above. The screen will appear as below:



Press the Exit key to return to the Supervisor Menu ↑ Press the General Options key ↓



Press the Exit Key to return to the General Options Menu ↑ Press the Change I.D.s & Passcodes Key ↓



This screen displays keys corresponding to each level of access. Levels of access exceeding that of the I.D. and Passcode used to access the supervisor menu will be inaccessible.

Pressing a key will allow I.D.s and Passcodes at that level of access to be changed.




The screen displays a list of all I.D.s and Passcodes at the selected access level. By pressing an I.D. or Passcode data entry field a keypad will appear and enable the user to change that I.D. or Passcode.

Pressing the Reset Level I.D./Passcodes key will present the user with a warning and ask them to confirm their decision. Upon doing so all of the I.D.s and Passcodes at this level of access will be reset to their default values.

Press the Exit key to return to the General Options Menu.

# CYCLE LOG STORE MENU


The Cycle Log Store Menu can be accessed by pressing the Cycle Log Store Menu key from the Auxiliary Menus screen and inputting a valid I.D. and passcode combination. The screen will appear as below:



Press the Internal/External Storage key to switch between accessing files stored on the system's internal memory or on files stored on a USB drive. Key will display Internal Storage symbol or External Storage symbol as appropriate.

Press the Exit key to return to the Auxiliary Menus screen.

Press the View Last Cycle Run or enter the desired cycle number to view the data log for that cycle.




Use the arrow keys to navigate the cycle log. The small arrows keys will move the data log by one line in the appropriate direction, the large Up and Down arrow keys will navigate to the top or bottom of the display respectively.

Pressing the Delete Cycle Log key will prompt the user to confirm their decision. On doing so the selected cycle log will be deleted.

Press the Exit key to return to the Auxiliary menu screen.

Once the printing operation has been cancelled or has completed the Cycle Log Store Menu will be displayed.

Press the P key to print the selected cycle log.



# BOILER CONTROL MENU

## Integral Boiler Control

The Boiler Control Menu can be accessed by pressing the Boiler Control Menu key from the Auxiliary Menus screen. When the boiler is not running, the screen will appear as below:

| Boiler Control Menu                             |        |                 |                      |                    |          |
|---|--------|-----------------|----------------------|--------------------|----------|
|   |        | CURRENT STATUS  |                      | AT LAST LOCKOUT    |          |
| BOILER IS RUNNING                               |        | NO              |                      | NO                 |          |
| WATER IS AVAILABLE TO FEED BOILER               |        | YES             |                      | YES                |          |
| WATER IS WITHIN THE CONTROLLED HIGH LEVEL BAND  |        | NO              |                      | NO                 |          |
| WATER LEVEL IS ABOVE THE MINIMUM LOW LEVEL      |        | NO              |                      | NO                 |          |
| WATER PUMP OVERLOAD HAS TRIPPED                 |        | YES             |                      | YES                |          |
| BOILER IS IN LOCKOUT MODE                       |        | YES             |                      | YES                |          |
| BOILER OVER PRESSURE DETECTED – PRESSURE SWITCH |        | NO              |                      | NO                 |          |
| BOILER OVER PRESSURE DETECTED - TRANSDUCER      |        | NO              |                      | NO                 |          |
| BOILER LOW WATER LEVEL DETECTED                 |        | NO              |                      | NO                 |          |
| BOILER OVER TEMPERATURE DETECTED                |        | NO              |                      | NO                 |          |
| BOILER MCB SHUNT TRIP HAS BEEN OPERATED         |        | YES             |                      | YES                |          |
| BOILER MCB HAS TRIPPED (BANK 1, 2, 3)           |        | NO NO NO        |                      | NO NO NO           |          |
| OVERFILL/WATER LEVEL SENSOR ERROR               |        | NO              |                      | NO                 |          |
| EMERGENCY SWITCH HAS BEEN OPERATED              |        | NO              |                      | NO                 |          |
| AIR SERVICE FAILED                              |        | NO              |                      | NO                 |          |
| FILLING WATCHDOG TIMEOUT                        |        | NO              |                      | NO                 |          |
| HEATING WATCHDOG TIMEOUT                        |        | NO              |                      | NO                 |          |
| HEATER CONTACTOR FAILED CLOSED (BANK 1, 2, 3)   |        | NO NO NO        |                      | NO NO NO           |          |
| HEATER CONTACTOR FAILED OPEN (BANK 1, 2, 3)     |        | NO NO NO        |                      | NO NO NO           |          |
| START BOILER                                    |        | EXIT            |                      |                    |          |
| CHAMBER PRESSURE                                | 0959mb | BOILER PRESSURE | 4.5bar               | JACKET TEMPERATURE | 135.3°C  |
| JACKET PRESSURE                                 | 3160mb |                 |                      | LAST CYCLE NUMBER  | 000101   |
| DOOR 1 LOADING SIDE                             |        | DATE            | Wednesday 09/05/2018 | TIME               | 10:26:28 |

The screen displays boiler information that can be used for fault diagnosis in the event of a boiler lockout.

Press the Start Boiler key


| Boiler Control Menu                             |        |                 |                      |                    |          |
|---|--------|-----------------|----------------------|--------------------|----------|
|   |        | CURRENT STATUS  |                      | AT LAST LOCKOUT    |          |
| BOILER IS RUNNING                               |        | YES             |                      | NO                 |          |
| WATER IS AVAILABLE TO FEED BOILER               |        | YES             |                      | YES                |          |
| WATER IS WITHIN THE CONTROLLED HIGH LEVEL BAND  |        | YES             |                      | NO                 |          |
| WATER LEVEL IS ABOVE THE MINIMUM LOW LEVEL      |        | YES             |                      | NO                 |          |
| WATER PUMP OVERLOAD HAS TRIPPED                 |        | NO              |                      | YES                |          |
| BOILER IS IN LOCKOUT MODE                       |        | NO              |                      | YES                |          |
| BOILER OVER PRESSURE DETECTED – PRESSURE SWITCH |        | NO              |                      | NO                 |          |
| BOILER OVER PRESSURE DETECTED - TRANSDUCER      |        | NO              |                      | NO                 |          |
| BOILER LOW WATER LEVEL DETECTED                 |        | NO              |                      | NO                 |          |
| BOILER OVER TEMPERATURE DETECTED                |        | NO              |                      | NO                 |          |
| BOILER MCB SHUNT TRIP HAS BEEN OPERATED         |        | NO              |                      | YES                |          |
| BOILER MCB HAS TRIPPED (BANK 1, 2, 3)           |        | NO NO NO        |                      | NO NO NO           |          |
| OVERFILL/WATER LEVEL SENSOR ERROR               |        | NO              |                      | NO                 |          |
| EMERGENCY SWITCH HAS BEEN OPERATED              |        | NO              |                      | NO                 |          |
| AIR SERVICE FAILED                              |        | NO              |                      | NO                 |          |
| FILLING WATCHDOG TIMEOUT                        |        | NO              |                      | NO                 |          |
| HEATING WATCHDOG TIMEOUT                        |        | NO              |                      | NO                 |          |
| HEATER CONTACTOR FAILED CLOSED (BANK 1, 2, 3)   |        | NO NO NO        |                      | NO NO NO           |          |
| HEATER CONTACTOR FAILED OPEN (BANK 1, 2, 3)     |        | NO NO NO        |                      | NO NO NO           |          |
| STOP BOILER                                     |        | EXIT            |                      |                    |          |
| CHAMBER PRESSURE                                | 0959mb | BOILER PRESSURE | 4.5bar               | JACKET TEMPERATURE | 135.3°C  |
| JACKET PRESSURE                                 | 3160mb |                 |                      | LAST CYCLE NUMBER  | 000101   |
| DOOR 1 LOADING SIDE                             |        | DATE            | Wednesday 09/05/2018 | TIME               | 10:26:28 |

Press the Exit key to return to the Auxiliary Menus screen.



## Printer

The printer is situated on the front control panel. The printer provides a record of process variables and times throughout the sterilizing cycle – the cycle log. In addition to the cycle log, the contents of the Supervisor, Validation and Configuration menus can be printed for archiving purposes.



### STATUS LED

The printer incorporates an LED indicator to report its condition. If there is a fault, the LED will flash in sequence. The fault can be identified by counting the number of flashes.

| LED INDICATION  | CONDITION              | SOLUTION           |
|-----------------|------------------------|--------------------|
| ●               | Printer On             | -                  |
| ●               | Printer Off            | -                  |
| ● ● ●           | Paper out or door open | Fit new paper      |
| ● ● ● ● ● ● ● ● | Thermal head too hot   | Allow head to cool |

Depressing the mode button will allow paper to be fed through the printer. When removing printout from the printer, pull the printout toward the tear bar and tear from one side to the other across the serrated edge.

### HEAD THERMAL LIMIT

After extensive printing the print head temperature may rise to an unusable level. The status LED will report when this occurs, and printing will be suspended until the head temperature returns to normal levels.

### PRINTER SELF-CHECK

A printer self-check can be performed to check the print quality of the printer.

- a) Turn the printer off.
- b) Push in and hold the mode button. Whilst still held in, push the on/off button to turn the printer on again.
- c) When the mode button is released, paper will start to be fed from the printer and characters to be printed. After 5 seconds release the paper feed button. The printer will continue to print until the self check is complete. Examine the characters for print quality and for correct operation of the print mechanism.

### REPLACING THE PRINTER PAPER

- a) The printer will automatically detect when the printer paper has run out, and report this using the Status LED. Use the Mode button to feed through the last centimeters of paper and fit a new roll.
- b) Pull the lever until the lid is released from its locked position. To avoid damage do not use excessive force.
- c) If the paper roll needs replacing, open the paper cup lid and remove the remaining paper. Reel off a few centimeters from a new roll of paper. Hold approximately 5cm of paper outside the device as you place the new roll into the reservoir. Close the lid by applying equal amount of pressure on each side ensuring the lid is in the locked position. Now tear the spare paper away.

## SERVICE REQUIREMENTS

### Spare Parts

The following spare parts are available from BMM Weston;

| Description                              | Part Number    | Qty   |
|--|----------------|-------|
| Printer "thermal" Paper                  | H15746         | 10    |
| Door Seal (Series 1)                     | H23270/H24398* | 2     |
| Door Seal (Series 2)                     | H12192/H21590* | 2     |
| Door Seal (Series 3)                     | H21441/H23776* | 2     |
| PTFE Silicon Grease                      | C7080          | 1 can |
| Air Filter                               | 60101          | 2     |
| O Ring - 6mm                             | 190623         | 10    |
| O Ring - 15mm                            | 190625         | 10    |
| Tri-Clamp Seal - 1"                      | H9189          | 10    |
| PTFE Tri-Clamp Seal - 1.1/2"             | H9190          | 10    |
| Seal Kit for steam inlet valves          | DS586          | 2     |
| Actuator Kit for steam inlet valve       | DS 575         | 2     |
| Seal Kit for Vacuum Valve                | DS 589         | 1     |
| Actuator Kits Vacuum Valve               | DS 576         | 1     |
| Seal Kit - 3/8" SMC Valve                | DS 767         | 2     |
| Spares Kit for Jacket steam Trap         | DS 510         | 2     |
| Push fit pneumatic elbow (1/4"BSP x 4mm) | C4598          | 5     |
| Fuse 1A                                  | H10010         | 2     |
| Fuse 2A                                  | H10011         | 2     |
| Fuse 3.15A                               | H10388         | 2     |
| Fuse 5A                                  | H10012         | 2     |

\*Contact BMM Weston Ltd. to determine appropriate door seal for your machine.

# SERVICE REQUIREMENTS

**FOR SPARE PARTS CONTACT BMM SERVICE DIVISION:**

BMM Service Helpdesk on 01795 597127  
or Main Switchboard on 01795 533441 then key 1 for Service Helpdesk  
BMM Service Fax No. 01795 538891  
BMM Service E-mail [service@bmmweston.co.uk](mailto:service@bmmweston.co.uk)  
BMM Weston Website [www.bmmweston.com](http://www.bmmweston.com)

Before contacting BMM Weston it would be useful to know

|              |                      |
|--------------|----------------------|
| Machine Type | <input type="text"/> |
| Serial No.   | <input type="text"/> |
| Model No.    | <input type="text"/> |

***THIS MANUAL WAS PRINTED AND PUBLISHED IN ENGLAND BY:-***

BMM WESTON LTD, WESTON WORKS, FAVERSHAM, KENT, ENGLAND

Head Office: Phone 01795 533441 Fax 01795 538891 Email: [mail@bmmweston.com](mailto:mail@bmmweston.com) Website: [www.bmmweston.com](http://www.bmmweston.com)

Service Dept: Phone 01795 597127 Fax 01795 538891 Email: [service@bmmweston.co.uk](mailto:service@bmmweston.co.uk)

**BMM WESTON**

Sterilizers